

Anti-His Polyclonal Antibody

Description

His-tag is widely used for cloning recombinant protein in different expression systems. Antibody against his-tag can be used as a common epitope to detect the expression of interested protein from recombinant cloning.

Product Name

Rabbit Anti-6xHis Polyclonal Antibody

Size: 100ug

Reactivity: human

Immunogen: 6x-His Peptide

Host: Rabbit

Applications

ELISA: 1: 1000- 1: 5,000.

Determine optimal working dilutions by titration test.

Storage: -20°C

Formulation

Antibodies are purified by protein A affinity chromatography

For Research Use Only