


Androgen Receptor (Phospho-Ser650) Antibody

Catalog No.	65441
Product Name	Androgen Receptor (Phospho-Ser650) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of androgen receptor only when phosphorylated at serine 650.
Reactivity	Human
Immunogen	Synthesized phosphopeptide derived from human androgen receptor around the phosphorylation site of serine 650 (T-T-S ^P -P-T).
Applications	IHC: 1/50 – 1/100
Swiss-Prot No.	P10275 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze / thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site.

Images

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Androgen Receptor (Phospho-Ser650) antibody.

Left: p-peptide -
Right: p-peptide +

For Research Use Only