



BCL-2(Ab-70) Antibody

Catalog No.	75381
Product Name	BCL-2(Ab-70) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total BCL-2 protein.
Reactivity	Human
Immunogen	Synthesized non-phosphopeptide derived from human BCL-2 around the phosphorylation site of serine 70 (R-T-S ^P -P-L).
Tested Applications	WB: 1/500 – 1/1000; IHC: 1/50 – 1/100; ELISA
Entrez Gene No.	596 (Human)
Swiss-Prot No.	P10415 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BCL-2 (Ab-70) antibody.

For Research Use Only