


BRCA1(Ab-1524) Antibody

Catalog No.	75460
Product Name	BRCA1(Ab-1524) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total BRCA1 protein.
Reactivity	Human, Mouse
Immunogen	Synthesized non-phosphopeptide derived from human BRCA1 around the phosphorylation site of serine 1524 (Y-P-S ^P -Q-E).
Tested Applications	Western Blot: 1/500 – 1/1000; IHC: 1/50 – 1/100; ELISA
Entrez Gene No.	672 (Human); 12189 (Mouse)
Swiss-Prot No.	P38398 (Human); O15129 (Human); Q1RMC1 (Human); Q3LRJ0 (Human); Q3LRJ6 (Human); Q6IN79 (Human); Q7KYU9 (Human); P48754 (Mouse)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using BRCA1(Ab-1524) antibody.

For Research Use Only