


ATM (Ab-1981) Antibody

Catalog No.	75468
Product Name	ATM (Ab-1981) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total ATM protein.
Reactivity	Human, Mouse
Immunogen	Synthesized non-phosphopeptide derived from human ATM around the phosphorylation site of serine 1981 (E-G-S ^P -Q-S).
Applications	IHC: 1/50 – 1/100
Swiss-Prot No.	Q13315 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze / thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using ATM (Ab-1981) antibody.

Left: peptide –
Right: peptide +

For Research Use Only