



AFX (Ab-197) Antibody

Catalog No. 75483
Product Name AFX (Ab-197) Antibody
Size 50ug, 100ug
Concentration 1 mg/ml
Specificity Detects endogenous levels of total AFX protein.
Reactivity Human, Mouse, Rat
Immunogen Synthesized non-phosphopeptide derived from human AFX around the phosphorylation site of serine 197 (A-A-S^P-M-D).
Applications Western Blot: 1/500 – 1/1000; IHC: 1/50 – 1/100
Entrez Gene No. 4303 (Human); 54601 (Mouse)
Swiss-Prot No. P98177 (Human); Q9WVH3 (Mouse); Q4KL34 (Mouse); Q3TZJ4 (Mouse)
Host Species Rabbit
Clonality Polyclonal
Isotype IgG
Storage -20°C for one year. Avoid freeze / thaw cycles.
Purification Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue, using AFX (Ab-197) antibody.

Peptide - +

For Research Use Only