


AMPK1 (Ab-174) Antibody

Catalog No.	75624
Product Name	AMPK1 (Ab-174) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total AMPK1.
Reactivity	Human, Mouse, Rat
Immunogen	Synthesized non-phosphopeptide derived from human AMPK1 around the phosphorylation site of threonine 174 (L-R-T ^P -S-C).
Applications	IHC: 1/50 – 1/100
Entrez Gene No.	5562 (Human); 105787 (Mouse); 65248 (Rat)
Swiss-Prot No.	Q1313 (Human); Q5EG47 (Mouse); P54645 (Rat)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze / thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images

Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using AMPK1 (Ab-174)antibody. Left: p-peptide – Right: p-peptide +

For Research Use Only