


FER (Ab-402) Antibody

Catalog No.	78100
Product Name	FER (Ab-402) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total FER protein.
Reactivity	Human, Mouse
Immunogen	Synthesized non-phosphopeptide derived from human FER around the phosphorylation site of tyrosine 402 (V-N-Y ^P -E-E).
Applications	WB: 1/500 – 1/1000
Swiss-Prot No.	P16591 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Buffer	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



For Research Use Only