

DATA SHEET

FOR RESEARCH USE ONLY

Catalog No.: C0011

Product Name: Human gastric adenocarcinoma (grade I) tissue array

Species: Human Rat Mouse Rabbit


Tissue Type: Adult Fetal Tumor Normal Disease Adjacent

Design Description:

- Each two cores come from one paraffin tissue block
- Adenocarcinoma, tubular adenocarcinoma, mucous adenocarcinoma. See chart below.
- Row : Column = 6 : 10
- Positive control: normal placenta (1); Negative control: blank paraffin (1)

Cases	Tissues	Tissue Blocks	Diameter	Thickness
29	1	29	1.5 mm	3 - 5 μ m
Related Products	CM005, C0011, C0027, C0032			

Array Layout:



Donor ID:

Donor Information:

Note: The clinical information provided above: for informational purposes only. Full Moon BioSystems, Inc.: not liable for any liability caused by the use of this information. If you need further assistance, please feel free to contact us.