

Array Scanner Recommendations

All scanners that are compatible with 75mm x 25mm (3in x 1in) microscope slides can be used to scan Full Moon BioSystems' antibody arrays.

Recommended Scanning Resolution: 20um or higher (10um, 5um, etc.)

Selected Compatible Systems

Manufacturer	Product Name	Required Accessory
Agilent Technologies	DNA Microarray Scanner	
Alpha Innotech	AlphaScan™ Microarray Scanner NovaRay™ Detection Platform	
Applied Precision	arrayWoRx® Biochip Reader	
Aurora Photonics	PortArray 5000™	
Biomedical Photometrics	The DNAscope™ HR	
Eppendorf North America	Eppendorf Silverquant® scanner system	
GE Healthcare, formerly Amersham Biosciences	Typhoon™ 9410, 9210, 8610	Microarray Slide Holder Kit
Genewave	AmpliReader™ 4600 Microarray Reader	
INNOPSYS	InnoScan® Microarray Scanner 700, 900, Autoloader	
Molecular Devices, formerly Axon Instruments	GenePix® Microarray Scanner 4000A, 4000B, 4100, 4200L	
PerkinsElmer, formerly GSI Lumonics	ProScanArray® HT Microarray Scanner ProScanArray® Microarray Scanner ScanArray® GX Microarray Scanner ScanArray® GX PLUS Microarray Scanner	
Tecan	LS Reloaded™ Versatile Scanner PowerScanner™ High Resolution microarray scanner	
Vidar Systems	Revolution™ 4200 Microarray Scanner	

Incompatible Systems

Manufacturer	Product Name
Affymetrix	GeneChip® Scanner HT Array Plate Scanner
Illumina	BeadStation 500GX Genetic Analysis System